Instruction for Use

Light Cured Dental Adhesive

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Light Cured Dental Adhesive

General Information

EsBond® is the 5th generation single component bonding agent designed to bond all classes of direct composite restorations to dentin, enamel, as well as composite, treated porcelain and set amalgam.

Composition

- Organic Matrix** : Pyromellitic dimethacrylate
2-Hydroxyethyl methacrylate
Acetone
- Additives** : Initiator, Catalyst, Stabilizer, Others

Indications

- All direct composite restorations
- Indirect restorative cementation: metal, porcelain and composite (For use with EsBond® Activator)
- Adhesive (bonded) amalgam restorations (For use with EsBond® Activator)
- Dual/self cure composite and core build-ups (For use with EsBond® Activator)
- Repairs of porcelain, composite, amalgam
- Composite to metal/set amalgam (Direct veneering)

Precautionary Information

- Check the expiration date and appearance condition of product and package.
- Observe the instruction for use after being well acquainted with the instruction for use.
- Do not use if the material has been cured.

Application

- Isolation** : Rubber dam is the preferred method of isolation.
- Tooth preparation and cavity preparation** : Rinse and dry. Remove minimal tooth structure. Prepare cavity with minimal tooth reduction. Use adequate pulp protection. Use a liner/base in areas of deep excavation.
- Etching** : Apply FineEtch® etchant to enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
Note. 1) In the case of excessive drying, it is difficult to form a hybrid layer due to collagen breakage, and the bonding strength is lowered.
2) Excessive acid corrosion increases the permeability due to dermal tubular fluid, resulting in post-operative hyper-sensitivity.
- Adhesive** : Immediately after blotting, apply 2-3 consecutive coats of adhesive (EsBond®) to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator.
Gently air thin for 5 seconds to evaporate solvents.
Note. 1) Please use EsBond® : EsBond® Activator (1:1) for self-cured or dual-cured.
2) When the solvent is evaporated after applying the adhesive, the dentin surface is applied continuously until it is glossy.

- Light cured** : Light cure for 20 seconds when the intensity of light is 1200 mW / cm².

	Intensity Light (mW/cm²)	Cure Time (sec.)
1	600	30
2	1200	20

- Restorative placement** : Place the restoration material according to the manufacturer's instructions.

Application method for Dual/self cure composite and core build-ups and indirect restorative cementation

- Isolation** : Rubber dam is the preferred method of isolation.
- Tooth preparation and cavity preparation** : Rinse and dry. Remove minimal tooth structure. Prepare cavity with minimal tooth reduction. Use adequate pulp protection. Use a liner/base in areas of deep excavation.
- Etching** : Apply FineEtch® etchant to enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- Adhesive** : Place mix equal number of drops of EsBond® and EsBond® Activator. Immediately after blotting, apply 2-3 consecutive coats of adhesive (EsBond® : EsBond® Activator (1:1)) to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator.
Gently air thin for 5 seconds to evaporate solvents.
Note. When the solvent is evaporated after applying the adhesive, the dentin surface is applied continuously until it is glossy.

IFU#101102(R 01/23)

Manufacturer

SPIDENT CO., LTD. www.spident.co.kr
203 & 312, Korea Industrial Complex, 722, Gojan-Dong, Namdong-Gu, Incheon, Korea 405-821
TEL. +82-32-821-0071 FAX. +82-32-821-0074

EU representative

AR EXPERTS B.V. (tradename : CE Medical)
Amerlandseweg 7, 3621 ZC, Breukelen, the Netherlands
TEL. +31 (0)85 007 32 20

- Light cured** : Light cure for 20 seconds when the intensity of light is 1200 mW / cm².

	Intensity Light (mW/cm²)	Cure Time (sec.)
1	600	30
2	1200	20

- Restorative placement** : Place the restoration material according to the manufacturer's instructions.

Caution

- EsBond® is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
- EsBond® is only used for oral hard tissue.
- EsBond® contains methacrylic monomer. Some patients can cause allergic reaction to monomer.
- Do not use EsBond® for patients allergic to or hypersensitive to methacrylic monomers.
- If EsBond® causes an allergic reaction or oversensitivity, discontinue its use immediately.
- Avoid contact of EsBond® with eyes, mucosal membrane, skin and clothing.
- If EsBond® comes in contact with the skin or eyes, clean with pure water several times, and then consult with a medical specialist if necessary.
- Use a rubber dam to protect the throat.
- DO NOT use EsBond® after the indicated expiration date on the syringe.
- If refrigerated storage is used, allow product to reach room temperature prior to use.
- Avoid the use of Phenolic substances or other phenolic containing materials in combination with dental adhesives because Phenolic substances, especially preparations containing eugenol or thymol or clove oil, interfere with the curing of dental adhesive materials.

Storage & Use

- Do not expose restorative materials to elevated temperatures, or to intense light.
- Store EsBond® at a temperature between 2 - 27°C.
- Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life.
- Shelf life at room temperature is 36 months.
- See outer package for expiration date.
- EsBond® adhesive should be capped immediately after use it minimize evaporation.
- Do not store materials in proximity to products containing Eugenol or thymol or clove oil.

Warranty

Spident® warrants this product will be free from defects in material and manufacture. Spident® MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Spident®'s sole obligation shall be repair or replacement of the Spident® 's product.

Limitation of Liability

Except where prohibited by law, Spident® will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.
Caution : U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

Disposal

It should be handled according to local laws and regulations.

Conformity with Comunaute Europeenne			
	Temperature limit		Batch code
	Caution		Date of Manufacture
	Keep dry		Used-by date
	Keep away from sunlight		Consult Instructions for Use
	Manufacturer		Medical Device
	Authorized representative in the European community		
	Federal law prohibits dispensing without a prescription		