


✓ Multi-shade build-up guide with EsCom® 250




CLASS I & II

※ Source: Spident R&D center internal data

DENTIN				ENAMEL
Posterior Dentin Shade				Posterior Enamel Shade
A1	A2	A3	A3.5	TW

Case 1

» posterior direct resin filling
#36, 37 (class I & II)



※ Prof. DuckSoo Kim(Kyunghee Univ.)

Opaque Shade	Dentin Shade		Enamel Shade
A02	A1	A2	TW



CLASS III

※ Source: Spident R&D center internal data

OPAQUE DENTIN		DENTIN				ENAMEL
Opaque Shade		Dentin Shade				Enamel Shade
A02	A03	A1	A2	A3	A3.5	TW

Case 2

» anterior teeth diastema | #11 #12




※ Daseul Jung(Jung, J.K Dental Clinic)

Opaque Shade	Dentin Shade		Enamel Shade
A02	A1	A2	TW

Case 3

» anterior teeth restoration | #11 (class III)



※ Daseul Jung(Jung, J.K Dental Clinic)

Opaque Shade	Dentin Shade		Enamel Shade
A02	A1	A2	TW



EsCom® 250

Light-Cured Restorative Nano-Hybrid Composite Resin

Your perfect choice,
The world's first easy-slide system,
is in your one hand now.

EZ SLIDE
Syringe



CE
1639



ECMS05R00/210823



SPIDENT CO., LTD.
203 & 312, Korea Industrial Complex, 722, Gojan-dong,
Namdong-Gu, Incheon, Korea 405-821
TEL. +82-32-821-0071 FAX. +82-32-821-0074



Certification Experts B.V. (tradename: CE Medical)
Amerlandseweg 7, 3621 ZC, Breukelen, the Netherlands
TEL. +31 (0)85 007 32 20






✓ Sufficient physical properties


>> Wear Resistance

Because nano fillers disperse uniformly, asperity rate is very low.
Therefore the asperity does not interrupt the opposing tooth moving.



※ Source: Spident R&D center internal data


>> Mechanical Properties



※ All features were measured according to the ISO standards by R&D center of SPIDENT CO., LTD.

>> Shape Retention


Paste shapes comparison between before and after 3, 5 min in 37 deg C




※ Source: Spident R&D center internal data

✓ Practical design for easy-slide system

>> No more wasting





>> EZ SLIDE Cap



✓ Formulation with balanced Nano Filler

>> Compared with EsCom® 250 – FE-SEM


- Natural Esthetics
- Excellent Physical Properties



※ Source: Spident R&D center internal data

✓ Color adaptability for natural teeth

>> Each shade is aligned in the order of color saturation. Four shades are placed and cured.



※ Source: Spident R&D center internal data

